

AEM/AEL ELECTROLYSER TEST STATION



- RSI make AEM/AEL ELECTROLYSER TEST BENCH is fully designed and integrated with software controlled (Made in India) to meet complex research needs of the researchers.
- AEM electrolyser test bench is focussed to provide modularity in design and also complies with international standards from every aspect.
- RSI make AEM electrolyser test bench is mainly used for research purposes. The researchers can use this test bench to perform various experiments and tests, such as measuring the current-voltage characteristics, assessing the durability and lifetime of the membrane and electrodes, and optimizing operating conditions for maximum hydrogen production.
- Send us your requirements and we will prepare a detailed offer for your turn-key project.



SYSTEM FUNCTIONS:

| SYSTEM SPECIFICATIONS | |
|--|---|
| Model No | RSI-AE <mark>M-AEL-3</mark> 00W |
| Туре | AEM |
| Valuma Of Electrolytic tank | (Anion Exc <mark>hange Me</mark> mbrane) 1000 mL |
| Volume Of Electrolytic tank Active Area Of The Cell | 3 X 3 cm Or 5x5 |
| Weight | 400 Kg Approx. |
| Weight | (*Subjected to the model) |
| Dimension (L x BX H) | 800 X 70 <mark>0 X700</mark> mm |
| | (*Subjected to the model) |
| Storage Method | Tedlar bag |
| Safety Features | Emergency stop |
| PUMP SPECI | FICATIONS |
| Electrolyte flow | 20 mL to 250 mL/min |
| Speed variation facility | |
| Set Flow/ Volume and time | |
| Forward and reverse control | 4 |
| Manual / Software Controlled | Customer Requirement |
| POWER S | SUPPLY |
| Max Voltage | 10V |
| Max Current | 30 A |
| Power | 300 W |
| Techniques | OCP,CC,CV,CC-CV |
| Voltage sense | |
| Manual / Software controlled | Customer Requirement |
| MASS FLOW MET | ER (OPTIONAL) |
| Flow rate | 0 - 500 mL / min |
| Software Monitor | ✓ |
| | |



SYSTEM FUNCTIONS:


| ADDITIONAL FUNCTIONS | |
|--|---|
| Conductivity Measurement | X |
| Automatic Electrolyte Filling | X X |
| Line Heating And Monitoring | (i) Power(W) |
| Tank Heating And Monitoring | Voltage Voltage(V) 1 1/1 1/1 Power Voltage(V) 2 1/1 1/1 Current (A) 8 1/1 1/1 Report (A) 8 1/1 1/1 Report (A) 8 1/1 1/1 |
| Cell Heating And Monitoring | -24 -8 -79 -235 -7.8 |
| Electrolyte Tank Level Sensing Maintaining | g &23 -7.6 -7.6 -7.5 -7.4 \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ |
| O2 Flow Measurement | -21.5 \\ \frac{2}{3} -7.1 \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ |
| H2 Flow Measurement | -6.9 -20.5 -6.8 -6.7 |
| Back Pressure Controlling | -20 -6.6 -19.5 -6.5 |
| 8 Channel AUX Voltage Monitoring 322621 PM 032731 PM 032 | |
| 8 Channel AUX Temperature Monitoring H2 Leakage Indicator H2 Senso | or Value (0-100%) H2 Side Pressure Transmitter (Bar) O2 Side Pressure Transmitter (Bar) |
| Gas Filtration (Single Stage) | 0.0131226 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 |
| Hydrogen Leakage Detector | |
| SYSTEM REQUIREMENT | |
| Processor | Intel Core i5 or i7 |
| Memory (RAM) | 16 GB or higher |
| Storage | 512 GB SSD or larger |
| Graphics | Dedicated graphics card with at least 4 GB VRAM (NVIDIA GeForce or AMD Radeon) |
| CPII Operating System | Windows 11 |




Software Features:



- Techniques:
 OCP,Charge
- Sequencing Of Techniques
- Loop sequences into multiple cycles
- Computer-controlled peristaltic pump (optional)
- Log all data into excel format
- Real time Graphical Representation
- Export specific graph data
- Capacity Monitor (Ah Vs Time Vs Wh)
- Voltage Monitor (Voltage Vs Time)
- Temperature Monitor(channel 1-8)
- Hydrogen Production Monitor (mL/min Vs Time)





THANK YOU!



Please consult with our technical team to meet all your customised requirements

RESEARCH SUPPORTERS INDIA

FOLLOW US



ISO 9001: 2015

MANUFACTURING AND R&D HQ
No.63/1, Club House Road,
Maruthi Nagar, Hosur- 635109 ,Tamil Nadu ,India
Mail id :sales@rsindia.net, Ph: +91 8148274261,